Manufacturer:
Stryker Trauma AG
 Bohnackerweg 1
 2545 Selzach
 Switzerland
 phone: +41 (0) 32 641 66 66
 fax: +41 (0) 32 641 66 60


English

Maintenance of Torque Limiters


To ensure accurate torque measurement over time, Torque Limiters require periodic testing.

Torque Limiters		Must be tested every
	702749 – 1.5 Nm Torque Limiter (3mm system)	100 surgeries or 1 year
	702750 – 4 Nm Torque Limiter (4mm system)	100 surgeries or 1 year
	702751 – 6 Nm Torque Limiter (5mm system)	100 surgeries or 1 year

- Use Torque Tester (ref 702755) by choosing the proper 3mm, 4mm or 5mm Adapter on either side.

- Assemble Torque Limiters 702750 and 702751 with their respective Screwdriver Bits (702753 & 702754) and T-Handles (702427 & 702430)

Ensure the red needle is at «zero», then introduce the Torque Limiter assembly into the adapter and turn it slowly clockwise until it clicks.

The red needle indicates the torque generated by the Torque Limiter. The measured torque must be within 10% of the nominal torque as shown in the table below. We recommend taking several measurements to ensure the testing procedure was done correctly.

Torque Limiters		Torque Specification
	702749 – 1.5 Nm Torque Limiter (3mm system)	1.35 to 1.65 Nm
	702750 – 4 Nm Torque Limiter (4mm system)	3.6 to 4.4 Nm
	702751 – 6 Nm Torque Limiter (5mm system)	5.4 to 6.6 Nm

Torque Limiters within specification can be used for surgery. Torque Limiters close to specification limit should be tested frequently to ensure remaining within specification. Torque Limiters outside specification must be discarded and replaced by new ones.

Important Notes

Torque Testers and its adapters are not designed to be steam sterilized. Testing procedure must be performed in a non-sterile environment. Torque Tester must be regularly recalibrated by an accredited calibration laboratory. Torque Tester must be used with care and not dropped.

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