


S2™ Tibial Nail

IM Nail Case Report #1



S2™ Tibial Nail

IM Nail Case Report #1

Equipment:

S2™ Tibial Nail

Patient's injury:

methaphyso-epiphyseal fracture of right tibia and distal $\frac{1}{3}$ fibula

From:

Dr. Gilbert Taglang

Chief of Emergency Service

Centre de Traumatologie et d'Orthopédie (CTO)

Strasbourg, France

Background Information:

The patient, a 24 year-old male, suffered a football injury resulting in a comminuted, distal $\frac{1}{3}$ tibia fracture, with an associated fracture of the fibula. He came into the emergency room at Centre de Traumatologie et d'Orthopédie (CTO) in Strasbourg after he had been temporarily immobilized with a cast.

Considering the location of the fracture and its displacement, along with the patient's age and activity level (young and active), it was decided to reduce the distal metaphyseal tibial fragment first, using one cannulated screw.

Then, intramedullary nailing was performed. Because of its very distal locking holes, the S2™ Tibial Nail was selected for fixation.

Preoperative X-Rays



Procedure

Step 1: Fracture Reduction

To perform the procedure, the patient was put on a traction table in the supine position. The fracture was reduced via closed means, by applying transcalcaneal traction. To restore the tibial pilon, a 5mm Asnis III screw, 55mm in length was placed, anatomically reducing the fracture.

Step 2: Nail Insertion

A para-tendinous incision was made, extending approximately 4cm distally from the patella. The entry point was located under fluoroscopic control. The cortex was penetrated with the curved awl and the ball tipped guide wire passed down the tibial shaft. Progressive reaming was performed, up to 12mm in diameter.

A 10mm diameter nail, 330mm in length, was then inserted.

Step 3: Locking

Proximal locking was performed using the Target Arm. Because of the very distal position of the fracture, only one 5mm screw was used for proximal locking.

Distal locking was performed in a free-hand manner. For better rotational stability of the distal fragment, the nail was first locked with an AP screw and then with one ML screw (the more distal one, at 5mm from the tip of the nail).

Comments

The proximal ML distal hole was not locked, because it was located exactly at the level of the fracture site.

The distal fragment was very well fixed and stabilized with one AP and one ML screw.

Postop X-Rays



Asnis Screw insertion



Proximal Locking



Distal AP locking

Postoperative Follow-Up, from 8 Weeks

Results

Immediate mobilisation was permitted, without weight bearing for the first 8 weeks. Then, up to 3 months, partial weight bearing was allowed, also depending on consolidation level noticed on the follow-up X-rays.

The patient was doing extremely well at the eight week follow-up appointment. The fracture was healing and he was able to perform sports activities, 6 months after the surgery.



Dr. Gilbert Taglang

Chief of Emergency Service
Centre de Traumatologie et d'Orthopédie (CTO)
Strasbourg, France

Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

EMS Equipment

Stryker Trauma GmbH
Prof.-Küntschers-Strasse 1-5
D-24232 Schönkirchen
Germany

www.trauma.stryker.com

Always refer to the package insert, product label and/or user instructions before using any Stryker product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with [™] designation are trademarks of Stryker.
Products referenced with [®] designation are registered trademarks of Stryker.

Copyright © 2005 Stryker
Printed in Germany